



OraStripdx™ is a rapid exam-room test for the early detection of hidden periodontal disease in dogs and cats without sedation or anesthesia.

It works by detecting thiols — the compounds produced by the anaerobic bacteria associated with periodontal disease.

The test strip is gently glided along the upper gums at the gingival margin. If thiols are present, the indicator pad changes color from white to yellow within 10 seconds, providing visible evidence of infection hidden beneath the gums.

Seeing is believing. When owners are educated about how OraStripdx works and they can see the color change, they're

more likely to understand the need for dental care and follow veterinary recommendations. That's why studied hospitals saw a 54% increase in dental procedures when using the thiol test in every wellness exam.

OraStripdx is backed by published peer reviewed clinical and business use case studies which support product use, including validation studies that the thiol-detection test disclosed active periodontal disease in its earliest stages. The studies are found on our website at OraStripdx.com.

Key Takeaways

- Easily administered by veterinary teams in the exam setting
- Rapid color change makes test easy-to-read
- Effective for both cats and dogs; recommended for use in every wellness exam
- Proven to detect periodontal disease in its earliest stages
- Proven to increase client compliance with veterinary dental recommendations
- Proven to increase practice revenues from dental procedures, radiology, and home care products
- A thiol test is part of WSAVA and AAHA dental guidelines.

In this study, the thiol-detection test provided a firm, objective basis for veterinarian recommendations of dental procedures, including professional dental cleanings and radiography, and home-care products. Routine use of the thiol-detection test therefore resulted in clear and unambiguous recommendations and high rates of compliance. This in turn resulted in substantially increased practice dental revenues*.

Read the full article
at [OraStripdx.com](https://www.orastripdx.com)



Product Specs:

40 tests per tube
Detailed use instructions included
Storage: room temperature
Shelf life: 24 months