

Designed  
for Accuracy,  
Reliability and  
Simplicity



The Cardell® Insight Diagnostic Monitor is designed to be compact, lightweight and durable with an intuitive interface to make it easy to use.

#### DESIGNED FOR VETERINARY CLINICS

Featuring veterinary-specific algorithms that provide quick, accurate and consistent readings, the Cardell® Insight Diagnostic Monitor is useful in all areas of a veterinary practice, including exam, treatment, surgery and recovery. It can monitor patients in the field or those being transferred within the clinic.

#### BP TECHNOLOGY

The Cardell® BP technology you already know, adapts to your patient and employs advanced motion suppression to help eliminate errors caused by shivering, panting or other movement. This advanced feature can result in more efficient, reliable and timely readings.

#### FITS ALMOST EVERYWHERE

The Cardell® Insight Diagnostic Monitor is designed to fit anywhere within your clinic—on a narrow shelf or table, mounted to the wall, or on a mobile stand for ultimate mobility.



#### YOU'RE SUPPORTED

Midmark is 100% committed to your training, and we offer the instruction you and your team need to become completely proficient on your new equipment. With our free clinical support hotline, online resources and the most comprehensive in-clinic training program, Midmark gives you access to the tools that can help you provide better care for your patients.

# CARDELL® INSIGHT DIAGNOSTIC MONITOR SPECIFICATIONS

Non-Invasive Blood Pressure (NIBP)	Cardell® Technology
Pulse Oximetry (SpO <sub>2</sub> )	Nellcor® Oximax®
PR	NIBP PR
Display	7-Inch LCD Resistive Touch Screen
Waveform Display	SpO <sub>2</sub> (8015)
PR Range	30-240
NIBP Modes	Manual, Auto, Stat
Tabular Trends	72 Hours
Parameter Alarms	3 Defaults (Cat, Dog, Horse)
Download to PC	USB Drive
Weight	4.3 LB
Battery Life	10 Hours
Voltage	100-240 V
Warranty	2 Years, Limited

## CARDELL® INSIGHT DIAGNOSTIC MONITORS

- 8014-001** Diagnostic Monitor with NIBP
- 8015-001** Diagnostic Monitor with NIBP, SpO<sub>2</sub>

## MONITOR MOUNTS

- 8008-005** Mobile Stand
- 9A456012** Wall Mount

Additional monitor mounting options are available.

## COMPANION ACCESSORIES

The accessories necessary for full monitoring capability are included with each monitor.

Additional monitoring accessories are available at [midmarkserviceparts.com](http://midmarkserviceparts.com).



## COMPREHENSIVE ANESTHESIA MONITORING TRAINING

Our RACE-approved training courses, provided by a VTS in anesthesia and analgesia, can teach you how to achieve the highest level of anesthesia monitoring excellence to improve patient care. To learn more, visit [midmarkacademy.com](http://midmarkacademy.com).

- VI-AMT** Anesthesia Monitoring 101
- VI-AMT2** Anesthesia Monitoring 201