



DENTALAIRE

Meeting and Exceeding the Needs of Veterinary Dentistry since 1988

www.Dentalaire.com

800-866-6881

**High-speed
Handpiece**

Highspeed Handpiece 400LED

A New Generation Air Powered LED

Features:

- Ceramic Bearings
- 3-port Coolant Spray
- Push Button Auto Chuck
- Light Weight: 124 grams
- Standard ISO 4 Hole Configuration
- Autoclavable: 135° C / 275° F
- 400,000 rpm
- Quiet Operation
- One Year Warranty



Four Hole
Connection



Ceramic Ball
Bearings



**Bright LED Lights
Built-in Generator
Air Powered**

Technical Specifications

Head Type	Standard
Rotation Speed	370,000-400,000 rpm
Drive Air Pressure	Minimum 35 psi
Chuck Type	Push Button
Burs	Friction Grip Short Shank ISO 1797 Type3 1.59-1.60
Optics	Direct-LED
Head Size W x H	.4389in x .4192in

No Fiber-Optic Electronics Needed.