



KVP RENEWABLES

BetterCollars™



You won't notice the difference, but the *planet* will.

Designed to address the plastic waste issue in the veterinary industry, KVP BetterCollars are made of biodegradable PLA derived from renewable food sources like corn, and are made to breakdown when in landfills.

Features

- 100% plant based and biodegradable
- Complete degradation when exposed to moisture and heat in 3-6 months
- Proven to biodegrade based on ASTM D5526 (simulating landfill condition)
- FDA (510k) Approved
- Full body coverage and cushioned neckline and rim with durable, biodegradable plastic snaps

Sizes

ITEM#	SIZE	NECK CIRC.	DEPTH
BC-S	Small	6"-9.5"	4.5"
BC-M	Medium	7"-12"	5"
BC-L	Large	10"-15"	7"
BC-XL	X-Large	12"-21"	9.25"
BC-XXL	XX-Large	15"-25"	12"



KVP's BetterCollars are made of PLA, which stands for Polylactic acid. Alternate versions of this plastic are also used in 3D Printing!



BetterCollars

Renewable E-Collars

KVP BetterCollars are our latest addition to the KVP Renewables initiative.

Made from BIO-RESIN Polylactic Acid (PLA) Technology and plant-based renewable sources such as corn, these collars look, feel, and perform exactly like standard e-collars and include a premium stitched cotton lining.

Frequently Asked Questions

- **How will these products begin to break down?**

The collars will only breakdown once exposed to heat and moisture over a prolonged period of time, 3 - 6 months.

- **Can I throw this away?**

Yes, these products can be thrown away in normal trash. At the landfill, these will disintegrate fully in 3 to 6 months.

- **Can these be recycled?**

Yes, you can fully recycle these.

- **Can I put these in my compost pile?**

No, unlike food scraps, these collars are only compostable in an industrial compost facility. Please consult your local recycling facility for composting options.

- **Will light affect these collars?**

No, exposure to light will not effect the product.

- **Can these collars be sterilized?**

Yes, but only using gas or gamma sterilization processes.

Environmental Tests

USDA Bio-Based (**Certified**)

Biodegradation (**ASTM 6400**)

Eco Toxicity (**ASTM 6868**)

Disintegration (**ASTM 5338**)

European Packaging (**EN13422**)



To see our research or request samples visit kvpvet.com/bettercollars

KVP BetterCollars and KVP Renewables are KVP International Inc. affiliated brands. For more information visit www.kvpvet.com.