

VMC
VICTOR MEDICAL COMPANY



BMVet ELITE SERIES

**MULTI-PARAMETER
TOUCH SCREEN
MONITORS**


**UNCOMPROMISING PERFORMANCE.
COMPREHENSIVE CONNECTIVITY.**

The all-new, best-in-class Bionet BMVet Elite Series line of multi-parameter touch screen monitors all feature advanced cardiac analysis – recognizing and alarming for the three most common lethal arrhythmias, and reliable and accurate SunTech non-invasive blood pressure technology. The BMVet Elite models are versatile and lightweight, making them ideal for dental and wet stations. Easily mount to a wall mount or rolling cart to save precious clinic space. All come standard with WiFi, built-in printer, and free BT-Link PRO and BT-Link Mobile telemetry and data management software.

**And Bionet gives you the industry's
best warranty – 4 years.**



**BM7Vet
ELITE**



**BM5Vet
ELITE**



**BM3Vet
ELITE**

VMC
VICTOR MEDICAL COMPANY

Victor Medical Company

800.888.8908 www.victormedical.com

Victor Medical Company

BMVet ELITE SERIES Multi-parameter Touch Screen Monitors



UNCOMPROMISING PERFORMANCE. COMPREHENSIVE CONNECTIVITY.

	BM7Vet ELITE 	BM5Vet ELITE 	BM3Vet ELITE 
Standard Parameters	ECG, SpO ₂ , NIBP (with SunTech technology), Respiration Rate, Temperature	ECG, SpO ₂ , NIBP (with SunTech technology), Respiration Rate, Temperature	ECG, SpO ₂ , NIBP (with SunTech technology), Respiration Rate, Temperature
Optional Parameters	Mainstream EtCO ₂ , Sidestream EtCO ₂ , Dual Gas (Anesthetic Gas + EtCO ₂), Invasive Blood Pressure (Dual Port), Advanced Arrhythmia Package (7 channel ECG w/ 5 lead cable)	Mainstream EtCO ₂ , Sidestream EtCO ₂ , Dual Gas (Anesthetic Gas + EtCO ₂), Invasive Blood Pressure (Dual Port), Advanced Arrhythmia Package (7 channel ECG w/ 5 lead cable)	Mainstream EtCO ₂ , Sidestream EtCO ₂ , Dual Gas (Anesthetic Gas + EtCO ₂)
Enhanced Connectivity	Includes BT-Link PRO Includes BT-Link Mobile WiFi and wired LAN ready	Includes BT-Link PRO Includes BT-Link Mobile WiFi and wired LAN ready	Includes BT-Link PRO Includes BT-Link Mobile WiFi and wired LAN ready
Screen	12" Color TFT LCD Touchscreen	10" Color TFT LCD Touchscreen	8" Color TFT LCD Touchscreen
Data Saving	168 Hours Trend Data Stored	168 Hours Trend Data Stored	168 Hours Trend Data Stored
Printer	Standard, Built-In	Standard, Built-In	Standard, Built-In
Data Management Software	BT Link (Windows) BT Link Mobile (iOS, Android) USB Data Port	BT Link (Windows) BT Link Mobile (iOS, Android) USB Data Port	BT Link (Windows) BT Link Mobile (iOS, Android) USB Data Port
Warranty	4 Years Parts and Labor 	4 Years Parts and Labor 	4 Years Parts and Labor 

Victor Medical Company

800.888.8908 www.victormedical.com