

# SonoMe *Vet*



## POINT OF CARE HANDHELD ULTRASOUND

### 8 models

that fit all your everyday ultrasound  
needs & budgets for starters & experts



#### Lightweight and 100% Portable

- Models weigh only 200g ~ 250g
- Embedded Wi-Fi module

#### Easy to Learn and Use

- On-screen operation with minimal buttons

#### Image Optimization

- Gain, TGC, Focus, Harmonic, Reject, Depth, Dynamic Range
- PRF & WF (in color mode)

#### Annotation

- On-screen with virtual keyboard

#### Cost Effective

- Less than a traditional ultrasound probe
- No activation. No subscription. No user fees.

#### Free Vet-Specific Mobile App

- iOS 8.0 or above / Vet patient admission

#### Measurement

- B, C, PDI mode: Length, Area, Circumference, Angle, Trace
- M mode: Heart Rate, Time, Distance
- PW mode: Velocity, Heart Rate(2), S/D, Depth

#### For All Patient Sizes

- Small (4.4lbs~), Medium/Large (~110.2lbs & more)









CONTACT BIONET TODAY 

to learn how SonoMe Vet can meet all your ultrasound needs – and budget

(714) 734-1760 | [Sales@BionetUS.com](mailto:Sales@BionetUS.com)

[www.sonome-vet.com](http://www.sonome-vet.com) | [www.bionetus.com](http://www.bionetus.com)

## MODELS & SPECIFICATIONS

Model	Budget	Image quality	Convex	Linear	Mode	Frequency	Element	Animal size
 <b>5CB</b> General usage	Minimum	Good	✓		BW	3.2~5Mhz	128 ele	S/M/L
 <b>10LB</b> General usage	Minimum	Good		✓	BW	7.5~10Mhz	128 ele	S/M
 <b>5C</b> Extended usage	Affordable	Good	✓		BW/COLOR	3.2~5Mhz	128 ele	S/M/L
 <b>10L</b> Extended usage	Affordable	Good		✓	BW/COLOR	7.5~10Mhz	128 ele	S/M
 <b>14L</b> Extended usage	Affordable	Good		✓	BW/COLOR	10~14Mhz	128 ele	S/M
 <b>H5C</b> Extended usage	Affordable	Better	✓		BW/COLOR	3.2~5Mhz	192 ele	S/M/L
 <b>H10L</b> Extended usage	Affordable	Better		✓	BW/COLOR	7.5~10Mhz	192 ele	S/M
 <b>H5C10L</b> Extended usage	Reasonable	Better	✓	✓	BW/COLOR	C: 3.2~5Mhz L: 7.5~10Mhz	192 ele	S/M/L

General usage: when you need BW at minimum budget

Extended usage: when you need BW and Color at affordable budget

Animal Size: S (Small 4.4lbs~); M/L (Medium/Large ~110.2lbs and more)

Specification subject to change without prior notice.

All SonoMe Vet probes come with 1 year warranty

## MOBILE APP



- **Vet-specific mobile app**  
Vet patient admission / Vet
- **Easy to use**  
One screen operation / Minimal buttons  
Short learning curve
- **Image Optimization, Measurement, Annotation**
- **Store and review**  
Directly on your device (JPG and MP4)
- **iOS 8.0 or above**

SOMV-BC-V1.01-0822

CONTACT BIONET TODAY 

to learn how SonoMe Vet can meet all your ultrasound needs – and budget

(714) 734-1760 | Sales@BionetUS.com

www.sonome-vet.com | www.bionetus.com