

CENTRE™

Topical Wound & Burn Dressing

CENTRE™

Topical Wound & Burn Dressing

KINETIC™ TECHNOLOGIES 961 Beasley Street, Suite 270 Lexington, Kentucky 40509
877.786.9882 (Toll Free) 859.258.9160 (Phone) 859.258.9188 (Fax)

FEATURES

Hyaluronate Sodium & PSGAG Sterile solution
5mg HA/mL-50mg/Vial
50mg PSGAG/mL-500mg/Vial

DESCRIPTION

Centre™ is a Hyaluronate Sodium (HA) Polysulfated Glycosaminoglycan (PSGAG) topical wound and burn dressing indicated for use on pressure injuries, venous stasis ulcers, first and second degree burns, superficial wounds, partial and full-thickness wounds, traumatic wounds, surgical wounds, autograft procedures and donor sites. For use on cats, dogs, and horses (not intended for food animal use).

DIRECTIONS

1. Cleanse wound/surgical site.
2. Withdraw the HA/PSGAG solution.
3. Dispense topically onto the wound site.
4. Cover wound with non-adherent dressing.
5. Change dressing as needed with repeated applications of Centre.™

PART NUMBER

Net Contents: 5 Vials 1004-07-01



www.kineticvet.com

