

High-speed Swivel LED Handpiece



DA400LEDMH

The High-speed Swivel LED Handpiece's macro-angle design and high color rendering index provide direct light on the face of the handpiece that is high in intensity and provides larger coverage.

No Fiber Optic Tubing

- Innovative micro-generator
- Generates its own light source
- Brighter and larger light than traditional fiber optic system

Swivel Rotation

- Smooth and quiet operation
- Helps reduce arm fatigue during each use
- Superior torque

Triple Water Spray

- Provides higher patient comfort
- Effectively cools the entire surface of the bur
- Helps avoid overheating

Durable Titanium Body

- Micro-textured and scratch-proof surface
- Solid titanium construction
- Built to last
- Reduces slippage

Compatibility

- Variety of coupling systems
- Highly adaptable to existing units

Technical Specifications

Head Size (W x H) mm	10.8 x 13.0
Power (WATT)	23
Drive Air Pressure (PSI)	40 - 60
Speed (RPM)	410,000
Water Spray	Triple Ports
Color Temperature (K)	5,500
Brightness (LUX)	25,500
Ø / max. length (MM)	1.6/25
Chuck Type	Push Button
Rotary Instruments	FG burs according to EN ISO 1797-1
Warranty	
Handpiece	1 year
Turbine	3 months