

E-Series I – Table Top Electric Dental Unit



Ultima Electric Dental Series

Ultima Dental and Bien-Air® have cooperated in developing the first electric driven dental machine specifically dedicated to small animal and equine dentistry. This breakthrough product provides the profession with a dedicated Ultima platform in conjunction with Bien-Air® Swiss electric handpiece technology. The compact system's electric handpiece delivers consistent torque differing significantly from air driven dental systems that can lose more than 50% of its speed when in contact with the surface to be cut.

Ultima Electric Dental Series features:

- Consistent handpiece torque
- No brushing techniques when contacting dental surface as the handpiece maintains speed regardless of pressure
- Rapid wrist motion is eliminated relating to hand, wrist and arm fatigue
- Complete system – new oil free mini-compressor used to cool handpiece and pressurize the water system
- No plumbing required with 2 - 1 liter water bottles

Ultima Electric Dental Series consists of:

- 1 - 40,000 RPM electric handpiece
- 1 - 1:1 straight nose cone for speeds of 100 to 40,000 RPMs
- 1 - 5:1 contra angle for speeds of 10,000 to 200,000 RPMs
- 1 Satelec® LED Piezo Scaler
- 1 - 3 way water, air and mist syringe
- Speed control to dial in your preferred speed

- **Precision** - dental cuts will be more precise as the cutting bur does not need to leave the cutting surface
- **No Handpiece Vibration**
- **Precise Speed** - using either the power dial or the variable speed foot control
- **Minimal Maintenance** - only need to oil and maintain the handpieces
- **Small Footprint**



Item # 54029

**CONSISTENT TORQUE -
NO BRUSHING TECHNIQUE
NEEDED.**

Item	Product Description	UOM
54029	E-Series I – Table Top Electric Dental Unit	Each