

# Olympic Elite Exam Lift Table™

**New – more efficient exam room table**  
**Saves time – Improves care – Patient friendly**

**Modern look for any exam room.**

**Corian® Top**  
Non-slip top. Patient friendly.  
Anti-drip edge.

**Fold-down Scale**  
Weigh during exams. Folds down.

**Free standing**  
Patient access from 4 sides.

**2 Accessory Sockets**  
One on each side.

41"  
↑  
**Lift Range**  
↓  
15"

**Electric lift** for smooth, quiet operation.

**Fully enclosed**  
Safer, cleaner, quieter.

**6 Electrical Outlets**



Scan to watch video demonstration

or visit us at [www.olyvet.com](http://www.olyvet.com)

**2 Foot Switches** on sides of table.



The Olympic Elite Exam Room Lift Table improves patient care and saves time on every examination. The electric lift provides perfect patient presentation for large and small animals, and the patient can be raised or lowered during the exam for better visualization and minor treatment.

In addition to height adjustment, the free-standing Elite Exam Lift Table provides complete access to the patient from all 4 sides for more thorough exams.

The Elite Exam Lift Table is also patient friendly and helps reduce anxiety. The warm feeling Corian® top is non-slip and non-reflective, so your patients can stand comfortably and calmly during an exam.

And with the Elite Exam Lift Table you no longer have to lift a heavy patient risking injury for your staff. It complies with ergonomic standards. For added safety and cleanliness, the electric lift column is fully enclosed.

With its Corian® top, electric lift and enclosed column, the Olympic Elite Exam Lift Table makes every exam room more efficient and welcoming.

## Accessories available

Choice of helpful Accessories mounted on Swing Arms, fit into sockets. No tools required.



The Olympic Elite Exam Lift is fully customizable for your procedural needs throughout the day. 2 sockets at the back of the table hold a sturdy 11"x20" Equipment Shelf with integrated hooks for stethoscope or clippers. 2 additional sockets can hold an Olympic IV Pole, LED Procedure Light, Swiveling Laptop/Monitor Tray or Instrument Tray. Because all accessories drop in to sockets, no tools are required to make changes.

# For your exam room – Adjustable height – 360° patient access



### Adjustable height

Better patient presentation for all patients large and small during any type of examination. Also, adjusts to comfortable standing or seated height for operator. Smooth, quiet electric lift has a capacity of 250 lbs. 2 foot switches, one on each side.



### Full patient access

The Olympic Elite Exam Lift table is free standing, not a built-in, so you have instant access to the patient from 4 sides for fast and thorough examinations. Saves time, improves patient care. Large Corian® surface, 54" x 22", height 37".



### Corian® – patient friendly

The Corian® top is warmer feeling and non-slip – patients stand comfortably and calmly during exams. Also, non-reflective so they aren't spooked by mirror-images and lights. No clawing at table. (See Corian® below.)



### No more lifting

With the Elite Exam Lift Table there is no more awkward lifting of heavy patients—and less likelihood of injuries for your staff. Large patients can walk on when the non-slip Corian® table is near floor level.



### Fully enclosed

The electric lift is fully enclosed in a telescoping column. Safer when other pets or children are in the exam room. And cleaner, with no exposed machinery to collect hair. Faster clean-up.



### Scale built-in – folds down

A precision scale is built in under the Corian® top. Scale has weight averaging, accuracy +/- 0.1 lb, lb/kg selection, large, easy-read numbers, and is swivel mounted for angle viewing. Lifetime accuracy guaranteed.

### Modern, attractive

The Olympic Elite Exam Lift Table makes any exam room look more attractive. The handsome Corian® top and sleek enclosed column give it an impressive modern look. Makes your hospital environment more pleasing to owners, pets and your staff.



### Corian® to match your décor

Corian® combines beauty with durability. It is the lifetime material that never ages from Dupont, not a thin laminate. Select from these standard color-patterns.



**Sahara**  
51280



**Sandstone**  
51281



**Platinum**  
51282

### Additional Features

**Casters** – The Elite Exam Lift Table rolls on 4 locking casters. Easily positioned in exam room or moved for cleaning floor. Flexible for future hospital changes.

**Rechargeable Battery** – Operates either connected to wall outlet or on internal fast-recharge battery.

**Power Strip** – Convenient power strip with four 120V outlets is mounted just below table top at rear. Handy for clippers and other devices.

## Order From Your Distributor

### Models with Corian® Table Top

**Olympic Elite Exam Lift Table - with Scale..... Model 51203**  
Includes rechargeable battery, Power Strip, sockets to hold up to 4 Olympic accessories and choice of Corian® Top.

**Olympic Elite Exam Lift Table - without Scale ..... Model 51204**  
Includes rechargeable battery, Power Strip, sockets to hold up to 4 Olympic accessories and choice of Corian® Top.

### Models with Stainless Steel Table Top

**Olympic Elite Exam Lift Table - with Scale..... Model 51205**  
Includes rechargeable battery, Power Strip and sockets to hold up to 4 Olympic accessories.

**Olympic Elite Exam Lift Table - without Scale ..... Model 51206**  
Includes rechargeable battery, Power Strip and sockets to hold up to 4 Olympic accessories.

### Accessories

**Equipment Shelf ..... Model 51613**

Drops into sockets next to scale

**Swiveling Instrument Tray, 11" x 13" ..... Model 51612**

**Ultra-Bright LED Procedure Light.....Model 50970**

**Advanced LED Procedure Light .....Model 50971**

**Swiveling Laptop/Monitor Tray ..... Model 51615**

**Drop-in IV Pole..... Model 51611**

**3-year Olympic Warranty. All parts and labor.**