

Olympic Fold-Up Gurney™

Saves space – Folding Gurney – Treatment Table



ROLLS



FOLDS



Compact storage
Only 8 inches deep.
Always ready for use.



Safety Lock
Automatically locks
when opened.



Scan to watch video demonstration

or visit us at www.olyvet.com

Features

Ideal for clinics with limited space – A transport gurney that folds out of the way and doubles as an exam-treatment table.

Folds and unfolds in seconds – To unfold, just spread the sides apart and an automatic safety lock holds it open. To close, release lock.

Rolls and steers easily – Large casters and steering handles at both ends. Turns on a dime.

Lightweight and strong – Weighs only 27 lbs but can carry 200 lbs. All welded construction.



Also a treatment table

When used with the Hard Top Stretcher, the Fold-Up becomes an exam-treatment table. Has firm, non-slip surface. Four locking casters.

Lift-off Stretchers

Safer. Easier. Faster.

The top is a stretcher held securely on 4 posts but easily lifts off for transferring or lifting patients. This lift-off stretcher helps reduce back injury and makes moving patients safer and faster. Eliminates awkward handling of compromised patients.



2 Stretchers available



Soft Stretcher

Very light-weight. Plastic fabric is medical grade, anti-microbial Herculite®. Wipes clean. X-ray transparent. Strong aluminum frame. Fully welded. Feet under corners of frame provide knuckle clearance.

Hard Stretcher

Converts gurney into treatment table. Molded of reinforced polymer. Non-slip surface. Wipes clean. Hard surface helps immobilize injured patients. Fluids drain to port with integrated bucket hook.

Order From Your Distributor

Olympic Fold-Up Gurney **Model 50604**
Complete with Hard and Soft Stretchers. Unfolded 52" x 21" x 36" tall.

Olympic Fold-Up Gurney **Model 50602**
With Soft Stretcher only.

3-year Olympic Warranty. All parts and labor.