

Olympic Ultra-Lite Gurney™

Storage space and lift-off stretchers



**EXAM
TREATMENT**



TRANSPORT



Scan to learn more

Lift-off Stretchers

Safer. Easier. Faster.

The top is a stretcher held securely on 4 posts but easily lifts off for transferring or lifting patients. This lift-off stretcher helps reduce back injury and makes moving patients safer and faster. Eliminates awkward handling of compromised patients.



Features

Versatile – 3-way use Here is an easy-roll gurney that can speed up patient handling, and be used for exams and treatments. Plus, it's got extra room for storage.

Ultra lightweight gurney – Because this gurney weighs only 31 lbs—just a fraction of what other gurneys weigh—it rolls effortlessly on its 4 swiveling casters.

Steering handles at both ends – So one or two people can steer and make tight turns in crowded spaces, narrow hallways.

Have supplies – will travel – Also with this gurney supplies can travel with the patient, so whatever you need is always with you.

Larger patients – The Ultra-Lite Gurney is 52 inches long and can carry patients up to 250 lbs. Heavy-duty construction – all welded.

Exam-treatment table – When used with the Hard Stretcher top, the Ultra-Lite becomes a treatment table with a reinforced surface that is firm and non-slip.

Extra storage – The generous storage space can hold up to 100 lbs of supplies and equipment (in addition to patient load).

MRI compatible model – A non-ferrous model of the Ultra-Lite is available for use in MRI environments.



**2 Stretchers
available**



Soft Stretcher

Very light-weight. Plastic fabric is medical grade, anti-microbial Herculite®. Wipes clean. X-ray transparent. Strong aluminum frame. Fully welded. Feet under corners of frame provide knuckle clearance.

Hard Stretcher

Converts gurney into treatment table. Molded of reinforced polymer. Non-slip surface. Wipes clean. Hard surface helps immobilize injured patients. Fluids drain to port with integrated bucket hook.

Order From Your Distributor

Olympic Ultra-Lite Gurney Model 50603
Complete with Hard and Soft Stretchers.

Olympic Ultra-Lite Gurney Model 50601
With Soft Stretcher only.

MRI Conditional Olympic Ultra-Lite Gurney Model 50690
With Hard and Soft Stretchers.

3-year Olympic Warranty. All parts and labor.