



#### VETERINARY INFUSION PUMP

- \* DIMENSIONS 5.7" X 4.7" X 3.9"
- \* WEIGHT 2.2 LBS
- \* PROGRAMMED ADMIN SETS: 10
- \* INFUSION HISTORY: 1,300 RECORDS
- \* RECALIBRATION: DONE WHILE IN THE HOSPITAL IN LESS THAN 10 MINUTES.
- \* WARRANTY: 24 MONTHS\*

\*WARRANTY FOR 24 MONTHS IS VALID FOR HOSPITALS THAT USE THE PRACTIVET IV ADMINISTRATION SETS.

#### VETERINARY SYRINGE PUMP

- \* LIGHTWEIGHT AND PORTABLE
- \* INCLUDES A CARRYING HANDLE
- \* STACKABLE - UP TO SIX
- \* STORES INFUSION HISTORY
- \* 5 INFUSION MODES: RATE, TIME, WEIGHT, INTERMITTENT MODE AND TOTAL IV ANESTHESIA MODE.
- \* FITS EITHER A VERTICAL OR HORIZONTAL POLE.
- \* APPROX 6 HOUR BATTERY LIFE

RECOGNIZES 5-60 ml MONOJECT, BD AND TERUMO SYRINGES

**PURCHASE EITHER A PRACTIVET INFUSION OR SYRINGE PUMP AT \$585 PLUS 4 CASES OF PRACTIVET IV CONSUMABLES (MIX AND MATCH) AND EARN A \$100 CREDIT\* ON THE PUMP FROM VICTOR MEDICAL SUPPLY.**

**\*CREDIT APPLIED IMMEDIATELY TO INVOICE.**



### STANDARD SET OPTIONS

VICTOR #	PRACTIVET #	Standard Administration Sets
1254004	AH7021	103" 10 Drop Admin Set w/ Inj Site and Rotating Male Luer
1254018	AH7117	110" 10 Drop Admin Set w/ 2 Inj Sites and Rotating Male Luer
1254005	AH7048	103" 15 Drop Admin Set w/ Inj Site and Rotating Male Luer
1254007	AH7054	78" 15 Drop Admin Set w/ Inj Site and Rotating Male Luer
1254006	AH7047	70" 60 Drop Admin Set w/ Inj Site and Rotating Male Luer
1254003	AH7019	103" 10 Drop Admin Set w/ Inj Site and Rotating Male Luer outside of cage
1254135	AH7053	32" (81cm) 10 Drop Admin Set w/ Rotating Luer, Non-DEHP Tubing

VICTOR #	PRACTIVET #	Standard Extension Sets
1254001	AH7043	30" M/F Ext Set w/ Inj Site and Rotating Male Luer
1254002	AH7207	42" M/F Ext Set w/ 2 Inj Sites and Rotating Male Luer
1254114	AH2102	50" Male/Male Extension Set
1254035	AH7039	60" Smallbore Ext Set 1.8ml Priming Volume w/ Slide Clamp, & RML
1254140	AH7061	36" Smallbore Ext Set .3ml Priming Vol w/ Clamp RML
1254139	AH7271	60" (152 cm) .42ml Priming Vol Smallbore Ext Set w/Clamp, Rotating Luer

VICTOR #	PRACTIVET #	Standard Catheter Connectors
1254130	AH7027	7" Smallbore Ext Set w/ T-Port Luer Slip
1254011	AH7107	7" Smallbore Ext Set w/ T-Port Luer Lock
1254009	AH7038	Male Luer Lock Inj Site (Male Adapter)
1254010	AH7090	Male Luer Slip Inj Site (Male Adapter)
1254017	AH7151	7" Regularbore Ext Set w/T-Port Luer Slip

VICTOR #	PRACTIVET #	BAXTER EQUIVALENT: Standard Administration Sets
1254029	AH8002	103" 10 Drop Admin Set Baxter pump Slide Clamp, Inj Y Site, & RML
1254030	AH8003	103" 10 Drop Admin Set Baxter pump Slide Clamp, 2 Inj Y Sites, & RML

VICTOR #	PRACTIVET #	BURETTE IV SET
1254041	AH7101	105" 60 Drop 150ml Burette IV Set

### NEEDLELESS SET OPTIONS

VICTOR #	PRACTIVET #	Needlefree Administration Sets
1254012	AH7181	80" 15 Drop Admin Set w/ NeedleFree Port, Roller Clamp, and Rotating Male Luer
1254015	AH7195	103" 15 Drop Admin Set w/ NeedleFree Port & RML
1254138	AH7171	100" 15 Drop Admin Set w/ Roller Clamp, 2 NeedleFree Ports, & Rotating Male Luer

VICTOR #	PRACTIVET #	BAXTER EQUIVALENT: Needlefree Administration Sets
1254027	AH8000	103" 10 Drop Admin Set, Baxter pump Slide Clamp, NeedleFree Port, & RML
1254029	AH8001	103" 10 Drop Admin Set, Baxter pump Slide Clamp, 2 NeedleFree Ports, & RML
1254141	AH8004	150" 10 Drop Admin Set Baxter pump Slide Clamp Inj site, MC port & RML

VICTOR #	PRACTIVET #	Needlefree and Specialty Extension Sets
1254137	AH7170	36" Extension Set w/ NeedleFree Port
1254025	AH7223	42" Male/Female Ext Set w/ 2 NeedleFree Ports, 2 Clamps, RML
1254032	AH7162	60" Smallbore Light Resistant Ext Set w/ Clamp & Luer Lock

VICTOR #	PRACTIVET #	Needlefree Catheter Connectors
1254013	AH7002	OmniSet 7" Smallbore NeedleFree Ext Set w/ T-Port Luer Lock
1254012	AH7001	OmniSet 7" Smallbore NeedleFree Ext Set w/ T-Port Luer Slip
1254015	AH3300	NeedleFree Male Adapter (Luer Lock)
1254014	AH7004	7" Smallbore Bifuse NeedleFree Ext Set, Clamps, Luer Lock
1254034	AH7145	6" Smallbore Trifuse Ext Set w/ 3 MicroClaves, 3 Clamps & Rotating Male Luer

VICTOR #	PRACTIVET #	Needlefree Bag Spikes and Vial Adapters
1254131	AH10	NeedleFree Bag Spike
1254116	AH50	NeedleFree Multi-Dose Vial Adapter
1254023	AH62	NeedleFree Connector with Vial Spike for 5cc Vial
1254132	AH70	Universal Vented NeedleFree Vial Access Device

VICTOR #	PRACTIVET #	Other Needlefree Products
1254134	AH7073	Smart Y Connector For Needle-Free Centesis (50 per Case)

**PLEASE CONTACT YOUR DISTRIBUTOR OR PRACTIVET  
PRODUCT SPECIALIST FOR MORE DETAILS**