


VSS® MILLENNIUM

25/30 KHz Ultrasonic Scaler

VSS2000

The Backbone to Your Veterinary Office

- Versatile and dependable
- Clearly marked intuitive controls
- Automatically adjusts frequency to power any 25KHz or 30KHz insert
- Autoclavable handpiece sheaths improve infection control
- On-demand foot-controlled Turbo Power Boost for heavy calculus
- Broad range of power easily accommodates ALL patients
- Manual purge and priming feature



- A. Convenience and control with the Foot Pedal.
- B. Preinstalled in-line water filter (VSS2040); one of two included with unit.
- C. The male quick-connect coupler at the end of the water hose allows you to easily connect the unit to your water source.
- D. Includes two autoclavable handpiece sheaths, made of high-temperature resin to last and keep asepsis (VSS2010)
- E. Includes two stainless steel inserts: (1) #10 30K Insert (VSS4068) & (1) Universal Perio 30K Insert (VSS4069).

Specifications

- **Power:** Input: 100-240 VAC, 1.0-2.0 A, 50/60 Hz; Output: 24 VDC, 3.0 A, 2.5mm, Center Positive
- **Weight:** 0.73 Kg / 1.60 lbs.
- **Size:** 1 ½" H x 5 ½" W x 7 ½" D
- **Weight of Transformer:** 0.6 lbs. (270g)
- **Length of Handpiece Cable:** 7' (215cm)
- **Length of Foot Pedal Cable:** 7' (215cm)
- **Length of Water Hose:** 7'6" (230cm)
- **Length of Power Cord:** 110V: 12'6" (380cm)
- **Handpiece Operating Frequency:** 25/30KHz
- **Protection Against Electric Shock:** Class 1, Type B applied part
- **Protection Against Ingress of Liquids:**
 - Foot Switch & Scaler: IPX1 (drip proof);
 - Power Supply: IPX0 (Ordinary)
- **Mode of Operation of Equipment:** Intermittent; 10 minutes on, 5 minutes off
- **Operating Conditions:** 15 - 30°C, 20 - 70% RH (non-condensing)
- **Transport and Storage Conditions:** -17°- 40°C, 20 - 70% RH (non-condensing)

ITEM #	PRODUCT
VSS2000	VSS® Millennium Scaler Control Unit with attached hard-wired handpiece, removable, autoclavable sheath (VSS2010), foot pedal, water line with male quick-connect coupler, and 2 Inserts.



Veterinary Services and Supplies, LLC

For more information, contact Veterinary Services and Supplies, LLC at **848-482-7201** or E-mail veterinaryservicesandsupplies@gmail.com