




VETCORDER™ PRO


EVERY PATIENT | EVERY TIME



Every Patient. Every Time.

Improving the quality of patient care has never been easier. VetCorder™ Pro is an affordable, handheld multi-parameter monitor capable of monitoring ECG, SpO₂ and internal temperature. With the available free app for Android and iOS, the VetCorder Pro can be mirrored via Bluetooth to a phone or tablet giving the veterinary team real-time, increased visibility of patient status from a distance. VetCorder Pro excels in many under-monitored areas of the clinic such as radiology, recovery and transport. With a 12-hour battery life, the pocket size VetCorder Pro is a must have for any clinic.





By providing enhanced monitoring capabilities, Vetcorder Pro improves patient safety and instills confidence in staff and animal owners. The pocket-sized touchscreen instrument does the job of multiple large monitors often exclusive to surgical suites. Affordable and economical at just \$995 (MSRP), Vetcorder Pro can help both new and established veterinary practices looking to optimize monitoring protocols and patient safety.

Monitoring Functions:

- Affordable patient monitoring in the palm of your hand
- Monitors SpO2, ECG and Internal Temperature
- Bluetooth connection with downloadable free app for Android and iOS

Patient care has never been easier!

Product Features:

- At a 4.5" x 2.5" x .75" size, it fits in the palm of your hand (and easily in your pocket)
- Weighs under 4 oz. – less than a typical smartphone!
- Bluetooth-enabled to allow continuous remote monitoring with a smartphone or tablet up to 30 feet away
- Battery longevity for up to 12 hours of continuous monitoring and can hold its charge for up to 3 months on standby



Uses:

Vetcorder Pro can be used on a wide variety of animal species, including dogs, cats, other companion animals, livestock, non-human primates, reptiles, birds and exotics.

It is ideal for:



Admission




Exam



Induction



Surgery



Dental



ICU



Transport



Recovery

Specifications Chart:

SpO2	
Range	70% - 100%
Accuracy	80-100% +/- 2%, 70-79 +/- 3%
HR/PR	
Range	30 to 350 bpm
Heart Rate	+/- 2 bpm or +/- 2%
Connectivity	
Wireless	Bluetooth 4.0
Battery	
ECG Recorder	Over 1,000 times
Bluetooth On	up to 12 hours
Battery Type	Rechargeable lithium-polymer 560 mAh
Vetcorder Pro Part Numbers	
MAI Animal Health	48561
MFG Part Number	7100-1750

Each System Includes:

- One Vetcorder Pro
- SpO2 lingual sensor
- ECG cable with two sets of clips
- USB charging cable with AC adapter
- Quick-start guide
- Warranty
- In-service training
- Internal temperature probe and extension

Using Vetcorder Pro from admission to discharge helps veterinary practices achieve a new level of patient security and peace-of-mind. The worst message to deliver to a client is the loss of a loved companion, especially when it could have been avoided. Vetcorder Pro keeps staff informed and connected to vital information to help ensure patients are not compromised.

Bluetooth connectivity via mobile device allows veterinarians, technicians and staff the flexibility to monitor patient vitals while executing other tasks or working with other patients in the practice or hospital.

Vetcorder Pro is widely used in:


- Veterinary practices and hospitals
- Exotic pet centers
- Emergency animal facilities
- Laboratories
- Zoos
- Universities
- Production animal environments

World Class Customer Support:

Every Vetcorder Pro system is backed by Sentier's dedicated service and support team. Purchase includes free in-service training with one of our technical experts and access to technical support via phone or web site. The Sentier team is always one call away to assist your practice or facility with questions, uses or technical support.

Bottom Line Patient Security:

Vetcorder Pro pays for itself in a short amount of time via minor monitoring fees. It is a rare occurrence that a pet parent would not want their loved family member both cared for and monitored from exam to recovery. It is the right thing to do for patient care.



The Case for Vetcorder Pro in Recovery:

According to a recent study, more than half of the deaths tied to veterinary surgery occurred within three hours of recovering from anesthesia, which is used in a multitude of minor procedures.¹

In a 2016 study, 57 percent of dogs and 35.5 percent of cats that underwent procedures using anesthesia experienced an adverse event, including: arousal or breakthrough pain (14.9%), hypoventilation (13.5%), hypotension (10.3%), arrhythmias (5.8%), hyperthermia/hypothermia (5.0%), airway complications (4.8%), recovery excitation (4.6%), aspiration risk (4.5%), desaturation (2.8%), hypertension (1.7%) and other ailments (3.7%).²

The study found that using monitoring tools like Vetcorder Pro can be easily integrated into veterinary practices to provide valuable evidence and data for veterinarians, veterinary technicians and anesthetists to maintain optimal health for their patients during procedures and in recovery.

¹Brodgelt, DC et al (2008) "The Risk of Death: The Confidential Enquiry into Postoperative Small Animal Fatalities." *Veterinary Anesthesia and Analgesia* 35, 365-373.

²McMillan, M et al (2016) "Adverse event surveillance in small animal anesthesia: an intervention-based, voluntary reporting audit." *Veterinary Anesthesia and Analgesia* 43(2): 128-35.



"We believe in the Vetcorder Pro. We're pet parents!"

This black Labrador retriever is being monitored with the Vetcorder Pro while recovering from a spay procedure. She is owned and loved by a Sentier employee. Wouldn't you want the best for your pets too?

Why Vetcorder?

Industry professionals say it best.

“Vetcorder does something no other monitor can do for us. It shows the status of our patients in recovery. Every other solution we’ve tried has fallen short of Vetcorder’s versatility and value. **With Vetcorder, our clinic can deliver a higher quality of patient care.**”

-Taylor Stewart, RVT

“After a week of using it, my team of **8 CVTs asked me to buy more units since they were fighting over the one I bought.** We now have a total of 3 and my other machines are gathering dust. The Sentier team has also been great about answering questions that we have. My techs and I highly recommend the Vetcorder.”

-Dr. Ilan Waskow, Mile High Veterinary Hospital

“Continuous monitoring is especially important when treating with drugs that affect the heart. Thanks to the Vetcorder, **I have peace of mind during transport and imaging.**”

-Shalee Ratliff, RVT

“We just can’t get enough of this device! **Use it in the ER, out in the field, post operatively or even during the intraoperative period. This device is incredibly versatile.** Did we mention its Bluetooth feature?”

-Veterinary Anesthesia Nerds “Top 10 Favorite Things” 2016, 2017 and 2018!

“We really enjoy using Vetcorder. It’s small size, ease of connectivity (even to our smartphones) and long battery life make it **an essential monitor in our practice before, during and in recovery on our anesthetic cases.**”

-Jan Bellows, DVM - President, All Pets Dental, Weston, Florida



Do the right thing for patients and clients.

Leverage the flexible, convenient power of Vetcorder Pro for patient monitoring from admission to discharge, ensuring animal safety and well-being. With Vetcorder Pro, it's easy to monitor.

Every patient. Every time!

To learn more, call your preferred distributor and set up a free two-week trial today or contact us at:
1-844-VETCORDER (1-844-838-2673)
info@sentierconnect.com
www.sentierconnect.com

©2020 MAI Animal Health. All rights reserved.



VETCORDER™ PRO
EVERY PATIENT | EVERY TIME

 **SENTIER**

 **MAI ANIMAL HEALTH**