



RICA Surgical Products



Weck® Horizon® Ligation Clip System

Precise Easy-Loading Titanium Clips

Weck's Horizon System provides veterinary surgeons all of the benefits and time-savings associated with the Hemoclip brand of ligation clips, PLUS these additional benefits:


- Easier to load clip, since the design of the Horizon cartridge creates a near frictionless one-hand load.
- Adhesive backing on the cartridge allows for placement in a convenient location.
- Another size clip for accurate and secure ligation of blood vessels of all sizes.

The Horizon System is the premium Weck titanium clip ligation platform, providing veterinary surgeons a strong, secure clip coupled with an efficient, convenient easy-load cartridge system.



The Horizon Clip

The pre-formed chevron shape of the Horizon clip engulfs tissue with precise tip-to-tip distal closure. The heart-shaped wire is designed to give each clip a firm grip on vessels, while the triangulated cross-section of the clip leg maximizes surface-to-surface contact between clip and jaw, practically eliminating clip fallout.



The Horizon Cartridge

The design of the Horizon cartridge creates a near frictionless one-hand load simplifying the loading of the clip into the applicator. Horizon cartridges are color-coded to match applicators and the adhesive backing allows convenient placement of the cartridge to any sterile area.

Save Time - Save Money

The ligation of blood vessels with ligation clips proceeds faster than the time to clamp and ligate vessels with conventional suture material. This results in a significantly shorter surgical time for the patient and the veterinarian.

Veterinarians report ligation clips are of great value for vessel ligation in the following types of surgical procedures:

- Feline and Canine Ovariohysterectomy
- Feline, Canine and Equine Castration
- Surgical Removal of Cutaneous and Subcutaneous Neoplasms
- Mastectomy
- Splenectomy, Nephrectomy, Adrenalectomy, Thyroidectomy
- Limb Amputation
- Tail Amputation, Toe Amputation
- Ocular Enucleation
- Intestinal Resection
- Patent Duct Arteriosus (PDA)
- Canine Vasectomy
- Liver Biopsy
- Persistent Right Aortic Arch
- Avian and Exotic Animal Surgical Procedures
- Equine Colic Surgery
- Plus many, many more...

See reverse side for the actual sizes and packaging of the sterile, titanium Horizon clips and current pricing.

Weck® Horizon® Ligation Clip System 2023 Veterinary Price List

LIST #	HORIZON CLIP SIZE	CLIPS/ CARTRIDGE	CARTRIDGES/ BOX	CLIPS/ BOX	VETERINARY PRICE PER BOX	VETERINARY PRICE PER CLIP
--------	-------------------	---------------------	--------------------	---------------	-----------------------------	------------------------------

Horizon Ligation Clip System – Titanium Clips

005200	Microclip™ (white)	6	30	180	\$217.00	\$1.21
001200	Small (yellow)	6	30	180	\$186.10	\$1.03
002200	Medium (blue)	6	30	180	\$219.80	\$1.22
003200	Medium/Large (green)	6	20	120	\$190.10	\$1.58
004200	Large (orange)	6	20	120	\$254.75	\$2.12

001204	Small (yellow)	24	25	600	\$566.90	\$0.94
002204	Medium (blue)	24	25	600	\$668.95	\$1.11
003204	Medium/Large (green)	24	15	360	\$517.30	\$1.44
004204	Large (orange)	24	15	360	\$693.30	\$1.93


LIST #	HORIZON APPLIERS	LENGTH	TYPE	PRICE
537061	Microclip Applier	6.0"	Curved (20°)	\$180.00
537081	Microclip Applier	8.0"	Curved	\$180.00
137061	Small Applier	6.0"	Curved	\$180.00
137081	Small Applier	8.0"	Curved	\$180.00
137085	Small Applier	8.0"	Angled (55°)	\$180.00
237061	Medium Applier	6.0"	Curved	\$180.00
237081	Medium Applier	8.0"	Curved	\$180.00
237085	Medium Applier	8.0"	Angled	\$180.00
337081	Medium/Large Applier	8.0"	Curved	\$180.00
337085	Medium/Large Applier	8.0"	Angled	\$180.00
337117	Medium/Large Applier	11.0"	Right Angled (70°)	\$180.00
437081	Large Applier	8.0"	Curved	\$180.00
437085	Large Applier	8.0"	Angled	\$180.00
437117	Large Applier	11.0"	Right Angled	\$180.00

Additional sizes of appliers are available. Contact RICA for further details.

DIMENSIONS OF WECK HORIZON CLIPS

HORIZON DESCRIPTION	WIDTH		HEIGHT		CLOSED CLIP LENGTH	
	INCHES	MILLIMETERS	INCHES	MILLIMETERS	INCHES	MILLIMETERS
Microclip Horizon (white)	0.075	1.91	0.075	1.91	0.088	2.24
Small Horizon (yellow)	0.078	1.98	0.113	2.87	0.143	3.63
Medium Horizon (blue)	0.119	3.02	0.184	4.67	0.232	5.89
Med-Large Horizon (green)	0.226	5.74	0.284	7.21	0.353	8.97
Large Horizon (orange)	0.307	7.80	0.396	10.07	0.488	12.38

NOTE: Width and height are taken from **INSIDE** the clip. For informational use only, all dimensions are **APPROXIMATE**.



To order, contact your authorized RICA Veterinary Distributor for current pricing, call RICA's office in Schiller Park, Illinois at 1-800-889-3218 or visit www.ricasurgical.com.



RICA Surgical Products

RICA Surgical Products, Inc.
9207 Ivanhoe Street
Schiller Park, IL 60176
(800) 889-3218

PRINTED IN USA
FORM # RSP-WHLCS-0923